

Establishing a Biorepository Collection

- Complete a [Biorepository Storage Request form](#) and submit to the Biorepository as described in the pdf.
- Your request for storage or access to specimens is reviewed, workflow and notifications are established prior to collection and storage.
- Questions may be directed to Torrey Gallagher at 603-650-7952

Placing an Order for the Biorepository

- LAB3175 Specimen to Pathology for tissue orders, complete “Biospecimen to Store” section
- LAB4043 Biorepository Request for non-tissue specimens (ie. Blood, stool, swab)
- Ensure the Study for which the specimen is being collected and stored is indicated either through the order dropdown selection or in the comments field of the order.

Guidelines for Biorepository Specimen Collection

- Patient name and MRN must be clearly written and attached directly to EACH specimen tube or container being submitted to biorepository.
- Labeling of the specimen tube or container is required to occur at the time of collection from the patient and that verification of the A# with the patient's wristband can be made.
- The requesting study is indicated on the Standard Laboratory requisition forms.

[Rothstein / WWC / 29424 / D17032 Study](#)

❖ Use a butterfly to draw this collection set

6mL lavender-topped (EDTA) blood tube

- Following venipuncture. invert the collection tube 10x
- Store at 4C (refrigerated) until courier pickup.
- Transport to the DHMC Main Campus Pathology LSR at Room Temperature
- LSR delivers to CGAT 4°C pass-through

10mL red-topped (Serum) blood tube

- Up to two (2) Red-topped Blood collection tube(s) are drawn
- Following venipuncture, invert the red-topped serum collection tubes 10x
- Allow the red-topped tubes to coagulate or clot at Room Temperature for 30 minutes.
- Transition coagulated serum tubes to 4C (refrigerated) until courier pickup.
- Transport to the DHMC Main Campus Pathology LSR at Room Temperature.
- LSR delivers to CGAT 4°C pass-through

10mL P800 blood collection tube **P800 Tubes Stocked at DHMC Heater Road Lab by WWC**

- Following venipuncture, immediately invert collection tubes and centrifuge for 20 minutes at 1300 rcf (3000 rpm).

- Maintain centrifuged tube in **upright** position following post-collection centrifugation, storage and transport to CGAT lab to avoid disrupting the blood separation.
- LSR delivers to CGAT 4°C pass-through

NOTE: Stool Collection Kit provided to Consented Patient by WWC, patient collects and returns to clinic, WWC RA deliver stool specimens to Heater Road for courier transport to main DH lab

10mL OMNI-Gut Stool Collection kit

- Stool Collector provided to patient by WWC study clinic
- Store at Room Temperature
- Specimen is provided to lab by patient; lab techs must label and document collection time
- Transport to the DHMC Main Campus Pathology LSR at Room Temperature
- LSR delivers to CGAT Biorepository Room Temperature

Gilbert-Diamond / WWC / 31064 / D19071 Study

- ❖ **MAIN CAMPUS COLLECTION SITE**
- ❖ **Blood Tube Study Kit provided to Phlebotomy by WWC RA at time of Blood Draw.**
- ❖ **Blood tubes labeled as usual at Blood draw**
- ❖ **Study blood collection transported by WWC RA to the ORO Lab for initial processing.**
- ❖ **Use a butterfly to draw this collection set**

Study Kit contains:

6mL lavender-topped (EDTA) blood tube
10mL red-topped (Serum) blood tube – up to two (2) tubes
10mL P800 blood collection tube

- Venipuncture completed using a butterfly needle
- Following venipuncture, tubes are inverted 10 times
- Study Kit tubes are provided to the WWC RA for transport to the DH ORO lab.

NOTE: Stool Collection Kit provided to Consented Patient by WWC, patient collects and returns to clinic, WWC RA deliver stool specimens to Heater Road for courier transport to main DH lab

10mL OMNI-Gut Stool Collection kit

- Stool Collector provided to patient by WWC study clinic
- Store at Room Temperature
- Specimen is provided to WWC Clinic by patient at Heater Road as described above

[Dragnev / ORIEN / 30172 / W16141 Study](#)

One 6mL lavender-topped (EDTA) blood tube

- following venipuncture, invert the collection tube 10x
- Store at 4C (refrigerated) until courier pickup.

- Transport to the DHMC Main Campus Pathology LSR at Room Temperature

[Tissue available in excess of clinical diagnosis](#)

Leach / T2M2 / 24281 / D13248 Study ENROLLMENT HAS ENDED

One 6mL lavender-topped (EDTA) blood tube

- following venipuncture, invert the collection tube 10x
- Store at 4C (refrigerated) until courier pickup.
- Transport to the DHMC Main Campus Pathology LSR at Room Temperature

One 10mL green-topped (Sodium Heparin) blood tube

- following venipuncture, invert the collection tube 10x
- Store at 4C (refrigerated) until courier pickup.
- Transport to the DHMC Main Campus Pathology LSR at Room Temperature

Tissue available in excess of clinical diagnosis

- [T2M2](#)
- [T2M2-GYN](#)
- T2M2-Lung
- T2M2-RCC
- [T2M2-Melanoma](#)